



Connect Your 3M™ True Definition to BioDent's Workflow

Save time and money when you send scans directly from your True Definition digital impression scanner. **In addition to all-ceramic restorations, BioDent accepts digital impressions for PFM and full-cast restorations.**

- Call 3M Digital service by dialing 800-634-2249. Select option 3, and then option 1.
- Request that they add BioDent to your account.
- Log out then log back into your 3M True Definition account.
- Select BioDent from the laboratory options.

If you have any questions, contact BioDent at 800-517-5250 and ask for digital impression technical support.

Looking for more information on our products?

A complete listing of restorative options and material descriptions, preparation and cementation guides and billing codes, visit biodentlaboratory.com/for-dentists

Save money on every digital submission
by choosing model free when possible.

SAVE \$20 On Your First Digital Case

Digital submissions are simple when you choose BioDent. Reference coupon code **GO_DIGITAL** on your Rx when you upload your first case and save \$20 on your first digital submission.

Coupon Code: GO_DIGITAL